

Product Datasheet

CD75 (B-Cell Marker); Clone LN1 (Ready-To-Use) Antibody (orb1712518)

Description	CD75 (B-Cell Marker); Clone LN1 (Ready-To-Use) Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Immunogen	Nuclei from pokeweed mitogen-stimulated peripheral blood lymphocytes.
Note	For research use only
Application notes	HeLa or Daudi Cells. Germinal center B-cells in a lymph node or tonsil. Use LN1's ability to stain surface of erythrocytes as an internal positive control.
Isotype	IgM, Kappa.
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.