

Product Datasheet

CD43 (T-Cell Marker); Clone 84-3C1 (Concentrate) Antibody (orb1712574)

Description	CD43 (T-Cell Marker); Clone 84-3C1 (Concentrate) Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Immunogen	Stimulated human leukocytes
Note	For research use only
Application notes	Paracortex in a tonsil or a reactive lymph node.
Isotype	IgG1, kappa
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.