

Product Datasheet

Smooth Muscle Myosin Heavy Chain (SM-MHC) (Leiomyosarcoma & Myoepithelial Cell Marker); Clone MYH11/923 (Concentrate) Antibody (orb1712782)

Description	Smooth Muscle Myosin Heavy Chain (SM-MHC) (Leiomyosarcoma & Myoepithelial Cell Marker); Clone MYH11/...
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Immunogen	Recombinant human MYH11 protein
Note	For research use only
Application notes	Uterus or normal breast.
Isotype	IgG1, kappa
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.