

Product Datasheet

CD45RA (Leukocyte Marker); Clone 158-4D3 (Concentrate) Antibody (orb1713249)

Description	CD45RA (Leukocyte Marker); Clone 158-4D3 (Concentrate) Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Immunogen	Stimulated human leukocytes
Note	For research use only
Application notes	Daudi, Raji, IM9 and human lymphocytes. Human lymph nodes and tonsils.
Isotype	IgG2a, kappa
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.