

Product Datasheet

Tau (pThr217) Antibody (HRP) (orb1713450)

Description	Mouse Anti-Human Tau (pThr217) Monoclonal IgG1 Kappa
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	HRP
Tested Applications	DOT, WB
Immunogen	Synthetic peptide of Wild-type Human Tau (pThr217) AA 210-224: SRTPSLP-pT-PPTREPK
Target	Tau (pThr217)
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% Sodium azide *Storage buffer may change when conjugated
Concentration	1mg/ml
Storage	-20°C
Note	For research use only
Application notes	A 1:1000 dilution of SMC-615 was sufficient for detection of Tau in 10 µg of Human, Mouse, Rat Tau Peptide (pThr217) by DB analysis using Goat Anti-Mouse IgG:HRP as the secondary antibody.
Isotype	IgG1 Kappa
Clonality	Monoclonal
Uniprot ID	P10636
Dilution Range	WB (1:1000); optimal dilutions for assays should be determined by the user.
Expiration Date	12 months from date of receipt.