

Product Datasheet

Alpha Synuclein A53T Mutant Pre-formed Fibrils (Type 1) (orb1714203)

Description

Human Recombinant A53T Mutant

formed Fibrils (Type 1)

Reactivity

Human

Tested

In vitro, In vivo, SDS-PAGE, WB

Applications
Target

Alpha Synuclein A53T Mutant

Preservatives

PBS pH 7.4

Concentration

Lot/batch specific. See included datasheet.

Storage

-80°C

Tag

No tag

Note

For research use only

Application notes

Certified > 95% pure using SDS-PAGE analysis. Low endotoxin 5 EU/mL @ 2mg/mL.

Protein Sequence

MDVFMKGLSK AKEGVVAAAE
KTKQGVAAEA GKTKEGVLYV
GSKTKEGVVH GVTTVAEKT
EQVTNVGGAV VTGVTAVAQK
TVEGAGSIAA ATGFVKKDQL
GKNEEGAPQE GILEDMPVDP
DNEAYEMPSE EGYQDYEPEA

Purity

> 95%

Source

Recombinant

MW

~14.46 kDa

Uniprot ID
P37840
Expiration Date

6 months from date of receipt.