


Product Datasheet

Rat Hey1(Hes-related family bHLH transcription factor with YRPW motif 1)ELISA Kit (orb1714653)

Description	Rat Hey1(Hes-related family bHLH transcription factor with YRPW motif 1)ELISA Kit
Reactivity	Rat
Range	31.25-2000pg/ml
Target	Hey1
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Rat
Hey1(Hes-
related
family
bHLH
transcr...