


Product Datasheet

Monkey VEGF(Vascular Endothelial Cell Growth Factor)ELISA Kit (orb1714670)

Description	Monkey VEGF(Vascular Endothelial Cell Growth Factor)ELISA Kit
Reactivity	Monkey
Range	31.25-2000pg/ml
Target	VEGF
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Monkey
VEGF(Vascular
Endothelial
Cell Gr...