


Product Datasheet

Mouse Oma1(Metalloendopeptidase OMA1, mitochondrial)ELISA Kit (orb1714677)

Description	Mouse Oma1(Metalloendopeptidase OMA1, mitochondrial)ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/ml
Target	Oma1
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse
Oma1(Metalloendopeptidase
OMA1, mi...