


Product Datasheet

Mouse 15-PGDH(Prostaglandin dehydrogenase 1)ELISA Kit (orb1714688)

Description	Mouse 15-PGDH(Prostaglandin dehydrogenase 1)ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	15-PGDH
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse 15-PGDH(Prostaglandin dehydrogenase 1)ELISA Kit