


Product Datasheet

Mouse GCLC(Glutamate-cysteine ligase catalytic subunit)ELISA Kit (orb1714691)

Description	Mouse GCLC(Glutamate-cysteine ligase catalytic subunit)ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/ml
Target	GCLC
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Mouse
GCLC(Glutamate-
cysteine ligase
cat...