

Product Datasheet

Mouse anti-SARS-CoV2(S-RBD)(Delta B.1.617.2)IgA ELISA Kit [Out of stock] (orb1714693)

Description	Mouse anti-SARS-CoV2(S-RBD)(Delta B.1.617.2)IgA ELISA Kit
Range	0.156-10ng/ml
Target	anti-SARS-CoV2 IgA
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.