


Product Datasheet

Mouse ACOD1/Irg1(aconitate decarboxylase 1)ELISA Kit (orb1714699)

Description	Mouse ACOD1/Irg1(aconitate decarboxylase 1)ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/ml
Target	ACOD1/Irg1
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Mouse
ACOD1/Irg1(aconitate
decarboxylase...