


Product Datasheet

Hamster anti- SARS-CoV-2(RBD)IgG ELISA Kit (orb1714701)

Description	Hamster anti- SARS-CoV-2(RBD)IgG ELISA Kit
Reactivity	Hamster
Range	1.563-100ng/ml
Target	anti- SARS-CoV-2 (RBD) IgG
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	0.938ng/ml
Expiration Date	Please enquire.



Hamster
anti-
SARS-
CoV-
2(RBD)IgG
ELISA K...