


Product Datasheet

Human KCNA3(Potassium voltage-gated channel subfamily A member 3)ELISA Kit (orb1714705)

Description	Human KCNA3(Potassium voltage-gated channel subfamily A member 3)ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	KCNA3
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
KCNA3(Potassium
voltage-gated
chan...