


Product Datasheet

Human EIF5A(Eukaryotic translation initiation factor 5A-1) ELISA Kit (orb1714708)

Description	Human EIF5A(Eukaryotic translation initiation factor 5A-1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	EIF5A
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
EIF5A(Eukaryotic
translation
initi...