


Product Datasheet

Human CLEC2D(C-type lectin domain family 2 member D)ELISA Kit (orb1714710)

Description	Human CLEC2D(C-type lectin domain family 2 member D)ELISA Kit
Reactivity	Human
Range	0.313-20ng/ml
Target	CLEC2D
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Human CLEC2D(C-type lectin domain family 2 member D)ELISA Kit