


Product Datasheet

Human Lactoferrin antibody IgG(LTF Ab IgG)Elisa kit (orb1714715)

Description	Human Lactoferrin antibody IgG(LTF Ab IgG)Elisa kit
Reactivity	Human
Range	1.563-100ng/ml
Target	LTF Ab IgG
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	0.938ng/ml
Expiration Date	Please enquire.



Human
Lactoferrin
antibody
IgG(LTF
Ab Ig...