


Product Datasheet

Human anti-SARS-CoV2(S-RBD)(Omicron,B.1.1.529)IgA ELISA Kit (orb1714726)

Description	Human anti-SARS-CoV2(S-RBD)(Omicron,B.1.1.529)IgA ELISA Kit
Reactivity	Human
Range	1.563-100ng/ml
Target	anti-SARS-CoV2(S-RBD) (Omicron,B.1.1.529)
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	0.938ng/ml
Expiration Date	Please enquire.



Human anti-SARS-CoV2(S-RBD)(Omicron, B.1...