


Product Datasheet

Monkey IFN- γ (Interferon Gamma)ELISA Kit (orb1714971)

Description	Monkey IFN- γ (Interferon Gamma)ELISA Kit
Reactivity	Monkey
Range	6.25-400pg/ml
Target	IFN- γ
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P63309
Sensitivity	3.75pg/ml
Expiration Date	Please enquire.



Monkey
IFN- γ (Interferon
Gamma)
ELISA Ki...