


Product Datasheet

Human SYNPO2(Synaptopodin-2)ELISA Kit (orb1714973)

Description	Human SYNPO2(Synaptopodin-2)ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	SYNPO2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9UMS6
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
SYNPO2(Synaptopodin-
2) ELISA Kit