


Product Datasheet

Rat SLIT3(Slit homolog 3 protein)ELISA Kit (orb1714977)

Description	Rat SLIT3(Slit homolog 3 protein)ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/ml
Target	SLIT3
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	O88280
Sensitivity	46.875pg/ml
Expiration Date	Please enquire.



Rat
SLIT3(Slit
homolog
3 protein)
ELISA ...