


Product Datasheet

Mouse SCN9A(Sodium channel protein)ELISA Kit (orb1714978)

Description	Mouse SCN9A(Sodium channel protein)ELISA Kit
Reactivity	Mouse
Range	0.313-20ng/ml
Target	SCN9A
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q62205
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Mouse
SCN9A(Sodium
channel
protein) ELIS...