


Product Datasheet

Human RAB3A(Ras-related protein Rab-3A)ELISA Kit (orb1714983)

Description	Human RAB3A(Ras-related protein Rab-3A)ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	RAB3A
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P20336
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
RAB3A(Ras-
related
protein
Rab-3A) ...