


Product Datasheet

Human HIOMT(Hydroxyindole O-methyltransferase)ELISA Kit (orb1715000)

Description	Human HIOMT(Hydroxyindole O-methyltransferase)ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	HIOMT
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P46597
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
HIOMT(Hydroxyindole
O-methyltransf...