


Product Datasheet

Human NUDT5(ADP-sugar pyrophosphatase)ELISA Kit (orb1715002)

Description	Human NUDT5(ADP-sugar pyrophosphatase)ELISA Kit
Reactivity	Human
Range	15.6-1000pg/ml
Target	NUDT5
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9UUK9
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
NUDT5(ADP-
sugar
pyrophosphatase)
E...