


Product Datasheet

Mouse NLRC5(Protein NLRC5)ELISA Kit (orb1715003)

Description	Mouse NLRC5(Protein NLRC5)ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/ml
Target	NLRC5
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	C3VPR6
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Mouse
NLRC5(Protein
NLRC5) ELISA
Kit