


Product Datasheet

Mouse Skp1(S-phase kinase-associated protein 1)ELISA Kit (orb1715005)

Description	Mouse Skp1(S-phase kinase-associated protein 1)ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/ml
Target	Skp1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9WTX5
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Mouse
Skp1(S-
phase
kinase-
associated
pro...