


Product Datasheet

Mouse OPA1(Optic Atrophy 1)ELISA Kit (orb1715009)

Description	Mouse OPA1(Optic Atrophy 1)ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/ml
Target	OPA1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P58281
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Mouse
OPA1(Optic
Atrophy 1)
ELISA Kit