


Product Datasheet

Human IL-22BP(Interleukin-22-binding protein)ELISA Kit (orb1715013)

Description	Human IL-22BP(Interleukin-22-binding protein)ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	IL-22BP
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q969J5
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human IL-22BP(Interleukin-22-binding pro...