


## Product Datasheet

### Guinea pig INS(Insulin) ELISA Kit (orb1715020)

|                        |                                   |
|------------------------|-----------------------------------|
| <b>Description</b>     | Guinea pig INS(Insulin) ELISA Kit |
| <b>Reactivity</b>      | Guinea pig                        |
| <b>Range</b>           | 7.813-500pg/ml □ 1ng=10uIU □      |
| <b>Target</b>          | INS                               |
| <b>Storage</b>         | 2-8 °C for 6 months               |
| <b>Note</b>            | For research use only             |
| <b>Uniprot ID</b>      | <b>P01329</b>                     |
| <b>Sensitivity</b>     | 4.688pg/ml                        |
| <b>Expiration Date</b> | Please enquire.                   |



Guinea pig  
INS(Insulin)  
ELISA Kit