


Product Datasheet

Human RAP1GAP(Rap1 GTPase-activating protein 1)ELISA Kit (orb1715029)

Description	Human RAP1GAP(Rap1 GTPase-activating protein 1)ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	RAP1GAP
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P47736
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
RAP1GAP(Rap1
GTPase-
activating
pro...