


Product Datasheet

Human CEL(Bile Salt-Activated Lipase)ELISA Kit (orb1715031)

Description	Human CEL(Bile Salt-Activated Lipase)ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	CEL/BAL
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P19835
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
CEL(Bile
Salt-
Activated
Lipase)
EL...