


Product Datasheet

Human MAMDC2(MAM domain-containing protein 2)ELISA Kit (orb1715044)

Description	Human MAMDC2(MAM domain-containing protein 2)ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	MAMDC2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q7Z304
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
MAMDC2(MAM
domain-
containing
prote...