


## Product Datasheet

### Human CD314(NKG2-D-activating NK receptor)ELISA Kit (orb1715049)

|                        |                                                     |
|------------------------|-----------------------------------------------------|
| <b>Description</b>     | Human CD314(NKG2-D-activating NK receptor)ELISA Kit |
| <b>Reactivity</b>      | Human                                               |
| <b>Range</b>           | 15.625-1000pg/ml                                    |
| <b>Target</b>          | CD314                                               |
| <b>Storage</b>         | 2-8 °C for 6 months                                 |
| <b>Note</b>            | For research use only                               |
| <b>Uniprot ID</b>      | <b>P26718</b>                                       |
| <b>Sensitivity</b>     | 9.375pg/ml                                          |
| <b>Expiration Date</b> | Please enquire.                                     |



Human  
CD314(NKG2-  
D-activating  
NK recepto...