


Product Datasheet

Human MCAT(Malonyl-CoA-acyl carrier protein transacylase, mitochondrial)ELISA Kit (orb1715051)

Description	Human MCAT(Malonyl-CoA-acyl carrier protein transacylase, mitochondrial)ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	MCAT
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	O43772
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
MCAT(Malonyl-CoA-acyl carrier prot...