


Product Datasheet

Human Paired Box Protein Pax-7(PAX7)ELISA Kit (orb1715055)

Description	Human Paired Box Protein Pax-7(PAX7)ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	PAX7
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P23759
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
Paired
Box
Protein
Pax-
7(PAX7)
ELI...