


Product Datasheet

Mouse anti-SARS-CoV2(S-RBD)(Delta B.1.617.2)IgG ELISA Kit (orb1715056)

Description	Mouse anti-SARS-CoV2(S-RBD)(Delta B.1.617.2)IgG ELISA Kit
Reactivity	Mouse
Range	1.563-100ng/ml
Target	anti-SARS-CoV2
Storage	2-8 °C for 6 months
Note	For research use only
Sensitivity	0.938ng/ml
Expiration Date	Please enquire.



Mouse anti-SARS-CoV2(S-RBD)(Delta B.1.61...