


Product Datasheet

Human cpeb1(Cytoplasmic polyadenylation element-binding protein 1)ELISA Kit (orb1715060)

Description	Human cpeb1(Cytoplasmic polyadenylation element-binding protein 1)ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	cpeb1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9BZB8
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
cpeb1(Cytoplasmic
polyadenylation ...