


Product Datasheet

Human BMPR2(Bone morphogenetic protein receptor type-2)ELISA Kit (orb1715067)

Description	Human BMPR2(Bone morphogenetic protein receptor type-2)ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	BMPR2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q13873
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
BMPR2(Bone
morphogenetic
protein r...