


Product Datasheet

Human AP1M2(AP-1 complex subunit mu-2)ELISA Kit (orb1715075)

Description	Human AP1M2(AP-1 complex subunit mu-2)ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	AP1M2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q9Y6Q5
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human AP1M2(AP-1 complex subunit mu-2) E...