


Product Datasheet

Human ALDH6A1(Malonate-semialdehyde dehydrogenase)ELISA Kit (orb1715087)

Description	Human ALDH6A1(Malonate-semialdehyde dehydrogenase)ELISA Kit
Reactivity	Human
Range	0.156-10ng/ml
Target	ALDH6A1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q02252
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Human
ALDH6A1(Malonate-
semialdehyde
dehy...