


Product Datasheet

Mouse COX4I1(Cytochrome c oxidase subunit 4 isoform 1, mitochondrial)ELISA Kit (orb1715092)

Description	Mouse COX4I1(Cytochrome c oxidase subunit 4 isoform 1, mitochondrial)ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/ml
Target	Cox4i1
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P19783
Sensitivity	0.094ng/ml
Expiration Date	Please enquire.



Mouse
COX4I1(Cytochrome
c oxidase subuni...