


Product Datasheet

Human SBDP 150(Alpha II SPECTRIN BREAKDOWN PRODUCT 150)ELISA Kit (orb1715093)

Description	Human SBDP 150(Alpha II SPECTRIN BREAKDOWN PRODUCT 150)ELISA Kit
Reactivity	Human
Range	0.313-20ng/ml
Target	SBDP 150
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q13813
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Human
SBDP
150(Alpha
II SPECTRIN
BREAKDO...