


Product Datasheet

Human SCN9A(Sodium channel protein type 9 subunit alpha)ELISA Kit (orb1715096)

Description	Human SCN9A(Sodium channel protein type 9 subunit alpha)ELISA Kit
Reactivity	Human
Range	0.313-20ng/ml
Target	SCN9A/Nav1.7
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q15858
Sensitivity	0.188ng/ml
Expiration Date	Please enquire.



Human
SCN9A(Sodium
channel
protein type ...