


Product Datasheet

Human CCT2(T-complex protein 1 subunit beta)ELISA Kit (orb1715108)

Description	Human CCT2(T-complex protein 1 subunit beta)ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	CCT2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P78371
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
CCT2(T-
complex
protein 1
subunit
b...