


## Product Datasheet

# Mouse Mydgf(Myeloid-derived growth factor)ELISA Kit (orb1715112)

|                        |                                                     |
|------------------------|-----------------------------------------------------|
| <b>Description</b>     | Mouse Mydgf(Myeloid-derived growth factor)ELISA Kit |
| <b>Reactivity</b>      | Mouse                                               |
| <b>Range</b>           | 15.625-1000pg/ml                                    |
| <b>Target</b>          | Mydgf                                               |
| <b>Storage</b>         | 2-8 °C for 6 months                                 |
| <b>Note</b>            | For research use only                               |
| <b>Uniprot ID</b>      | <b>Q9CPT4</b>                                       |
| <b>Sensitivity</b>     | 9.375pg/ml                                          |
| <b>Expiration Date</b> | Please enquire.                                     |



Mouse  
Mydgf(Myeloid-  
derived growth  
facto...