


Product Datasheet

Human EIF2B2(Translation initiation factor eIF-2B subunit beta)ELISA Kit (orb1715120)

Description	Human EIF2B2(Translation initiation factor eIF-2B subunit beta)ELISA Kit
Reactivity	Human
Range	31.25-2000pg/ml
Target	EIF2B2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	P49770
Sensitivity	18.75pg/ml
Expiration Date	Please enquire.



Human
EIF2B2(Translation
initiation fact...