


Product Datasheet

Human TUSC2(Tumor suppressor candidate 2)ELISA Kit (orb1715123)

Description	Human TUSC2(Tumor suppressor candidate 2)ELISA Kit
Reactivity	Human
Range	15.625-1000pg/ml
Target	TUSC2
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	O75896
Sensitivity	9.375pg/ml
Expiration Date	Please enquire.



Human
TUSC2(Tumor
suppressor
candidate 2...