


Product Datasheet

Human Anti-KLHL12(anti-Kelch-like protein 12 antibody)ELISA Kit (orb1715125)

Description	Human Anti-KLHL12(anti-Kelch-like protein 12 antibody)ELISA Kit
Reactivity	Human
Range	1.563-100ng/ml
Target	Anti-KLHL12
Storage	2-8 °C for 6 months
Note	For research use only
Uniprot ID	Q53G59
Sensitivity	0.938ng/ml
Expiration Date	Please enquire.



Human Anti-KLHL12(anti-Kelch-like protei...